



# Sample Handling



Solutions for sample processing  
in the clinical laboratory

## False bottom tubes - the optimal solution for laboratory automation



**NEW** Now available in 15 mm Ø and 13 mm Ø

- For direct analysis of decanted serum/plasma yields
- Straight wall and rounded base versions available to accommodate rack requirements

## Inserts



- Convenient for small volumes and direct analysis in analyzers
- Primary tube with identification is retained, preventing sample mix-up associated with secondary tubes

## StackPack Tubes



- Convenient, compact StackPack packaging
- 40% space savings compared to bulk packaging
- Available in a range of dimensions commonly required for secondary tubes


## Block racks - for efficient and economic sample storage




**NEW** Block rack with 13 mm Ø apertures

- Compact sample storage and organization
- Up to 9 colors provide coding options
- Available with or without handles
- Autoclavable and suitable for use in dishwashers
- Double block rack featuring 100 apertures available with 17 mm Ø (connected at lateral or front end)

## Anti-Evaporation Protection - hygienic and practical



**NEW** Universal cap: Can be used in two different positions to fit most tubes



**NEW** Universal anti-evaporation cap: Suitable for all standard blood collection systems or tubes from 13 to 16 mm Ø



■ Anti-evaporation caps for S-Monovette® 15 mm Ø



**NEW** Protective cover for block rack to protect open tubes (17 mm Ø) awaiting processing

## Tubes with round base in StackPack (without cap)

Volume	Order No.	Material	Length	Ø	Pcs./Inner box	Pcs./Outer box
5.0 ml	55.475	PS	75 mm	13.0 mm	500	2,000
5.0 ml	55.525	PP	75 mm	13.0 mm	500	2,000
5.0 ml	55.476	PS	75 mm	12.0 mm	500	2,000
5.0 ml	55.526	PP	75 mm	12.0 mm	500	4,000
8.0 ml	55.467	PS	100 mm	13.0 mm	500	2,000
11.5 ml	55.466	PS	100 mm	15.7 mm	250	1,000
11.5 ml	55.539.250	PP	100 mm	15.7 mm	250	1,000

Visit our 'Tube Finder' at [www.sarstedt.com](http://www.sarstedt.com) for a general overview of our total range of more than 250 tubes.

## Tubes with false bottom

- PP Tubes, HD-PE screw cap, packaged 100 pcs./bag, 1.000 pcs./box
- Printed version, features writing space, graduation and barcode line

Volume	Order No.	Length / Ø	Version	Tube base
3.5 ml	60.617.010	92 / 13 mm	without cap, without print	round
3.5 ml	62.617	92 / 13 mm	screw cap assembled, with print	round
3.5 ml	60.613.010	92 / 13 mm	without cap, without print	flat
3.5 ml	62.613	92 / 13 mm	screw cap assembled, with print	flat
5.0 ml	60.612.010	92 / 15.3 mm	without cap, without print	round
5.0 ml	62.612	92 / 15.3 mm	screw cap assembled, with print	round
5.0 ml	60.611.011	92 / 15.3 mm	without cap, without print, amber	flat
5.0 ml	60.611.010	92 / 15.3 mm	without cap, without print	flat
5.0 ml	62.611	92 / 15.3 mm	screw cap assembled, with print	flat

## Inserts

Order No.	Description	Compatible with	Packaging
73.1051	Insert 23.5 mm length	S-Monovette® Ø 13 mm	2,000 pcs./box
73.1052	Insert 23.5 mm length	S-Monovette® and tubes Ø 15 mm	2,000 pcs./box
73.1053	Insert 23.5 mm length	S-Monovette® and tubes Ø 16 mm	2,000 pcs./box
73.1052.001	Insert 32.0 mm length	S-Monovette® Ø 15 mm	2,000 pcs./box

## Anti-Evaporation Protection

Order No.	Description	Packaging
65.648	Anti-evaporation cap for S-Monovette® Ø 15 mm	5,000 pcs./box
65.1379	Universal cap for 12 to 16 mm Ø tubes (additional colours available on request)	10,000 pcs./box
65.647	Universal anti-evaporation cap for 13 to 16 mm Ø tubes	5,000 pcs./box
93.857	Protective cover for 17 mm Ø block rack, light grey	25 pcs./box

## Block Rack with 50 Apertures

- Our range includes up to 8 additional colours and multi-rack options. Brochures are available on request.

Aperture Ø	Version without handles, grey		Version with handles, grey	
	Order No.	Pieces / Box	Order No.	Pieces / Box
12.0 mm	93.854	20		
13.0 mm	93.853	16	93.855	12
17.0 mm	93.852	20	93.851	12