

# Thermo Scientific Finnpiquette C1

**The new Finnpiquette C1 is a portable, lightweight serological pipette holder. Designed for use with graduated and volumetric glass and plastic pipets from 1 ml to 100 ml, Finnpiquette C1 offers simple, efficient pipetting performance with maximum pipetting comfort.**

## Total control of pipetting

Finnpiquette C1 excels in offering total control of pipetting due to eight speed alternatives for both aspirating and dispensing. It also enables two different ways of “gravity only” dispensing. Additional control is gained from the 1 ml pipette support, that stabilizes pipetting with narrow 1 ml pipets. The Finnpiquette C1 also has a large illuminated LCD display, that shows battery charge, pump speeds for aspirating and dispensing, and low battery warning.

## Lithium-ion battery technology

Effective and comfortable to operate, Finnpiquette C1 features a heavy-duty lithium-ion battery that allows the pipettor to operate continuously approximately three times longer than regular units using nickel cadmium batteries. The lithium-ion battery can be recharged in just three hours. A universal power supply is delivered automatically with the unit.

## Compact storing between applications

Both a wall-mounted holder and a table-mounted stand are included in the package.



- Lithium-ion technology –provides up to 15 hours of continuous operating
- Lightweight ergonomic design
- Autoclavable silicon adapter
- 1 ml pipette support
- 2 ways of “gravity only” dispensing
- A universal power supply
- Fast recharging
- 2-year warranty with web registration

## Finnpiquette C1 Technical Specifications

Aspirate Speeds	8 options
Dispense Speeds	8 options (+gravity dispense)
Battery	Lithium-ion
Approx. Runtime Between Charges	15 hours of continuous use
Recharge time	3 hours
Weight	220 grams
Nosecone Pipet Holder	Autoclavable silicone rubber
Pipet Compatibility	All major brands of glass or plastic serological pipets, 1ml to 100 ml
Filter Assembly	0.45 micron standard
Filter Material	PTFE
Product specifications are subject to change without prior notice	

## Ordering information

Cat. no.	Description	Qty/Unit
4580800	Finnpiquette C1	1
<b>Filters</b>		
4580560	Filter, non-sterile	5/bag
4580570	Filter, sterile	5/bag
<b>Accessories and spare parts</b>		
22280	Pipette gripper housing	1
22278	Silicone pipette gripper	1
22300	Table mounted holder	1
22302	Wall Mounted holder	1